


12mm Annular Probes

Part Number	Description	Overall Dimensions			
		W	X	Y	Z
6531-52	Annular Electrodes 12 x 120mm Stainless Steel / Fortron	105	223	152	-
6531-52/PG	Pre fitted PG13.5 nut. Annular Electrodes 12 x 120mm Stainless Steel / Fortron	-	223	152	128
6532-52	Annular Electrodes 12 x 220mm Stainless Steel / Fortron	205	323	252	-
6532-52/PG	Pre fitted PG13.5 nut. Annular Electrodes 12 x 220mm Stainless Steel / Fortron	-	323	252	-
6530-52	Annular Electrodes 12 x 320mm Stainless Steel / Fortron	305	423	352	-
6530-52/PG	Pre fitted PG13.5 nut. Annular Electrodes 12 x 320mm Stainless Steel / Fortron	-	423	352	-
6535-52/PG	Pre fitted PG13.5 nut. Annular Electrodes 12 x 425mm Stainless Steel / Fortron	-	555	482	425
6534-52	Annular Electrodes 12 x 450mm Stainless Steel / Fortron	435	555	482	-
6534-52/PG	Pre fitted PG13.5 nut. Annular Electrodes 12 x 450mm Stainless Steel / Fortron	-	555	482	458



25mm Flush Probes


Part Number	Description	Approx Dimensions			
		W	X	Y	Z
6543-52	Probe Flush Electrodes 25 x 175mm Stainless Steel / Fortron	25	175	239	22
6541-52	Probe Flush Electrodes 25 x 75mm Stainless Steel / Fortron standard o-ring position	25	75	132	22
6542-52	Probe Flush Electrodes 25 x 75mm Stainless Steel / Fortron for use with some Braun saftey ports	25	75	132	46
6545-52	Probe Annular Electrodes 25 x 96mm Stainless Steel / Fortron standard o-ring position	25	96	154	22
6546-52	Probe Annular Electrodes 25 x 96mm Stainless Steel / Fortron for use with some Braun saftey ports	25	96	154	47

Futura Probes with certificates

Futura Probes with certificates - Add certificate Package Pt. No. 6520-64 for each probe purchased

FUTURA Probes

BIOTECH RANGE



Tri clamp probes

Part Number	Description	Overall Dimensions			
		W	X	Y	Z
6551-52	Probe Flush Electrodes 25 x 175mm Stainless Steel / Fortron with size 1 - 1.5" Tri-clamp flange, can be used with NovAseptic ports.	25	50.4 (1.984")	175	225
6552-52	Probe Flush Electrodes 25 x 70 mm Stainless Steel / Fortron with size 1 - 1.5" Tri-clamp flange, can be used with NovAseptic ports.	25	50.4 (1.984")	70	120
6553-52	Probe Flush Electrodes 25 x 70 mm Stainless Steel / Fortron with size 2" Tri-clamp flange, can be used with NovAseptic ports.	25	63.9 (2.516")	70	120

Futura Probes with certificates


Futura Probes with certificates - Add certificate Package Pt. No. 6520-64 for each probe purchased

FUTURA

Port Adaptors

BIOTECH RANGE

The port adaptors are designed to allow the 12mm straight Aber probes to be fitted into either Braun or Applikon 19mm (3/4") top plate, female port fitting. The New Brunswick version is similar but fits a male threaded port on either the top plate or side entry port. The o'rings provide a seal and the nut with nylon ferrules lock the probe in place.



Port Adaptors

Part No.	Description	Wetted O'ring 'a'	Sealing O'ring 'b'	Sealing Washer 'c'	Mating thread
1295-12	12mm probe to 19mm Braun Sartorius metric port	Silicone FDA 10.78 x 2.62	EPDM 26.7 x 1.78	-	26 x 1 mm PITCH
1296-12	12mm probe to 19mm Applikon imperial port	Silicone FDA 10.78 x 2.62	EPDM 26.7 x 1.78	-	3/4" BSP
1273-12	12mm probe to 19mm New Brunswick metric port	Silicone FDA 10.78 x 2.62	-	Silicone FDA 19 x 3.5 x 2	26 x 1 mm PITCH
1430-22	PG 13.5 Adaptor for 12 mm probes	Silicone FDA 10.78 x 2.62	-	-	PG 13.5