

FUTURA probes

12mm Annular probe

Product specification

Aber Instruments has a wide range of probes to suit all applications, the 7.5mm FUTURA pico probe for mini-bioreactors and similar smaller platforms, the 12mm diameter probes ideal for small bioreactors and 25mm diameter for larger vessels, typically used in manufacturing / production environments.

All probes are electro-polished and passivated. Options for USP class VI and FDA CFR21 177 certificates. By default a probe is shipped with a standard certificate of traceability. If you wish to have a full breakdown of all the materials used in the production of a probe then please order the 'Probe Certificate of Conformance' package for each probe. Contact ABER if you wish to see an example of the documentation you will receive.

For recommendations on suitable products for your application, please contact sales@aberinstruments.com



Annular Probe



Annular Probe
Pre fitted PG 13.5 Nut



FUTURA probes

12mm Annular probe

Product specification



Part Number	Description	Dimensions (mm)			
		W	X	Y	Z
6531-52	12 X 120mm Annular Probe	105	223	152	-
6532-52	12 X 220mm Annular Probe	205	323	252	-
6530-52	12 X 320mm Annular Probe	305	423	352	-
6534-52	12 X 450mm Annular Probe	435	555	482	-



Part Number	Description	Dimensions (mm)			
		W	X	Y	Z
6531-52/PG	12 X 120mm Annular Probe - Pre fitted PG 13.5 Nut	-	223	152	128
6532-52/PG	12 X 220mm Annular Probe - Pre fitted PG 13.5 Nut	-	323	252	228
6530-52/PG	12 X 320mm Annular Probe - Pre fitted PG 13.5 Nut	-	423	352	328
6534-52/PG	12 X 450mm Annular Probe - Pre fitted PG 13.5 Nut	-	555	482	458